

KEY ASSESSMENT AND MANAGEMENT OF

Mechanical Catheter Complications

FOR HEALTHCARE PROFESSIONALS ONLY The materials provided in this resource are for informational purposes only and are not intended as medical advice, or as a substitute for the medical advice of a physician/healthcare professional



CLICK HERE TO START





Obvious Pericatheter Leakage

This algorithm has been reviewed by:

Dr. Lim Chiao Yuen

Consultant Nephrologist Head of Nephrology Department Department of Renal Services Ministry of Health, Brunei Darussalam

Dr. Deane Cedraz

MD, MSc, Nephrologist Sr. Manager, Medical Affairs Renal Care (PD and HD) & Acute Therapies, Baxter APAC

External Catheter Tubing Damage



Baxter



Lower Than Expected Dialysate Drain Volumes

Outflow failure – Is it due to constipation?

FOR HEALTHCARE PROFESSIONALS ONLY The materials provided in this resource are for informational purposes only and are not intended as medical advice, or as a substitute for the medical advice of a physician/healthcare professional. AP-RC00-230003. 05/2023

Baxter

Lower Than Expected Dialysate Drain Volumes

+

bvious Pericatheter

External Cathe Tubing Damage PDACADEN

Lower T Expecte Drain Vo

Treatment

Oral laxatives

Review bowel habits with patient. Oral laxatives are needed to treat or prevent constipation which can impair the dialysate flow.

Baxter

vided in this resourc lical advice of a phys not intended Click here to 03. 05/2023 references

as medical advice,

Baxter



Lower Than Expected Dialysate Drain Volumes

NO Inflow good?

Inflow good?









Lower Than Expected Dialysate Drain Volumes bvious Pericatheter eakage

External Cathet Tubing Damage PDACADE

Lower T Expecte Drain Vo

Outflow failure – Is it du

CAUSES (P

Catheter malposition/ migration

Treatment options

- Radiologically guided manipulation
- Laparoscopically directed interventions
- Simultaneous catheter replacement

Baxter

Baxter

Lower Than Expected Dialysate Drain Volumes bvious Pericatheter eakage External Cathet Fubing Damage **PD**ACADEI



Baxter



Lower Than Expected Dialysate Drain Volumes bvious Pericatheter

External Cathet Tubing Damage **PD**ACADE











Lower Than Expected Dialysate Drain Volumes bvious Pericatheter eakage External Cathet Fubing Damage **PD**ACADEI

Lower T Expecte Drain Vo

Outflow failure – Is it due

CAUSES (pleas

Catheter kink



Treatment options

- Catheter manipulation via fluoroscopic guidewire
- Laparoscopic intervention
- Surgical manipulation

Baxter

Baxter



Lower Than Expected Dialysate Drain Volumes bvious Pericatheter eakage External Cathet Tubing Damage **PD**ACADEN

Lower T Expecte Drain Vo

Outflow failure – Is it du

CAUSES (ple

Intraluminal catheter occlusion (often by thrombus)



Treatment options

Flushing catheter with normal saline or heparinised saline

Baxter

Baxter

Lower Than Expected Dialysate Drain Volumes bvious Pericatheter eakage External Cathet Tubing Damage PDACADE

Lower T Expecte Drain Vo

Outflow failure – Is it (

CAUSES (p

Failure to dislodge intraluminal debris by brisk irrigation of catheter?

💼) Treatment options

 Thrombolytic therapy such as tissue plasminogen activator (tPA) for blood clot lysis

Baxter

Back

Baxter





Obvious Pericatheter Leakage

> * Leaks that occur ≤ 30 days post catheter implantation and start of PD

> > FOR HEALTHCARE PROFESSIONALS ONLY The materials provided in this resource are for informational purposes only and are not intended as medical advice, or as a substitute for the medical advice of a physician/healthcare professional. AP-RC00-230003. 05/2023

** Leaks that occur > 30 days post catheter

implantation and start of PD





Obvious Pericatheter Leakage

Early leak*

Dramatic leak from exit site and/or insertion incision?

Dramatic leak from exit site and/or insertion incision?

* Leaks that occur < 30 days post catheter implantation and start of PD | ** Leaks that occur > 30 days post catheter implantation and start of PD

FOR HEALTHCARE PROFESSIONALS ONLY The materials provided in this resource are for informational purposes only and are not intended as medical advice, or as a substitute for the medical advice of a physician/healthcare professional. AP-RC00-230003. 05/2023

Baxter

Lower Than Expected Dialysate Drain Volumes Obvious Pericatheter Leakage xternal Cathete ubing Damage PDACADEN



Click here to view references

Baxter



Obvious Pericatheter Leakage



Did the leak resolve with temporarily discontinuing dialysis for 1-3 weeks?

Did the leak resolve with temporarily discontinuing dialysis for 1-3 weeks?

* Leaks that occur < 30 days post catheter implantation and start of PD | ** Leaks that occur > 30 days post catheter implantation and start of PD



Lower Than Expected Dialysate Drain Volumes Obvious Pericatheter Leakage cternal Cathete Ibing Damage PDACADEN







Baxter

Lower Than Expected Dialysate Drain Volumes Obvious Pericatheter Leakage xternal Cathete ubing Damage PDACADEM



Baxter



Obvious Pericatheter Leakage

Late leak** Ultrasound or Computed Tomography (CT) scan

ULTRASOUND OR COMPUTED TOMOGRAPHY (CT) SCAN



Lower Than Expected Dialysate Drain Volumes Obvious Pericatheter Leakage kternal Cathete ubing Damage **PD**ACADEM







Lower Than Expected Dialysate Drain Volumes Obvious Pericatheter Leakage kternal Cathete ubing Damage PDACADEN



Baxter

FOR HEALTHCARE PROFESSIONALS ONLY The materials provided in this resource are for informational purposes only and are not intended as medical advice, or as a substitute for the medical advice of a physician/healthcare professional. AP-RC00-230003. 05/2023 Click here to view references



External Catheter Tubing Damage

Sufficient length of external tubing for repair?





External Catheter Tubing Damage



Repair external tubing, investigate peritonitis and treat as per clinic protocol. Peritonitis resolved?

Repair external tubing, investigate peritonitis and treat as per clinic protocol. Peritonitis resolved?



Lower Than Expected Dialysate Drain Volumes)bvious Pericatheter .eakage External Catheter Tubing Damage PDACADEN





FOR HEALTHCARE PROFESSIONALS ONLY The materials provided in this resource are for informational purposes only and are not intended as medical advice, or as a substitute for the medical advice of a physician/healthcare professional. AP-RC00-230003. 05/2023

Click here to view references

Lower Than Expected Dialysate Drain Volumes)bvious Pericatheter eakage External Catheter Tubing Damage **PD**ACADEI





FOR HEALTHCARE PROFESSIONALS ONLY The materials provided in this resource are for informational purposes only and are not intended as medical advice, or as a substitute for the medical advice of a physician/healthcare professional. AP-RC00-230003. 05/2023

Baxter



External Catheter Tubing Damage



With flow dysfunction / Satisfactory flow function and no concurrent peritonitis

FOR HEALTHCARE PROFESSIONALS ONLY The materials provided in this resource are for informational purposes only and are not intended as medical advice, or as a substitute for the medical advice of a physician/healthcare professional. AP-RC00-230003. 05/2023

Baxter

Lower Than Expected Dialysate Drain Volumes) bvious Pericatheter eakage External Catheter Tubing Damage **PD**ACADEN







Lower Than Expected Dialysate Drain Volumes)bvious Pericatheter eakage External Catheteı Tubing Damage **PD**ACADEN



Baxter

FOR HEALTHCARE PROFESSIONALS ONLY The materials provided in this resource are for informational purposes only and are not intended as medical advice, or as a substitute for the medical advice of a physician/healthcare professional. AP-RC00-230003. 05/2023

Click here to view references





References

- Elbokl, M., Momciu, B., Kishibe, T., Oliver, M. J., & Perl, J. (2020). Peritoneal dialysis access outcomes reported in randomized controlled trials: A systematic review. Peritoneal Dialysis International: Journal of the International Society for Peritoneal Dialysis, 089686082096689. https://doi.org/10.1177/0896860820966898
- Crabtree, J. H., Shrestha, B. M., Chow, K.-M., Figueiredo, A. E., Povlsen, J. V., Wilkie, M., Abdel-Aal, A., Cullis, B., Goh, B.-L., Briggs, V. R., Brown, E. A., & Dor, F. J. M. F. (2019). Creating and Maintaining Optimal Peritoneal Dialysis Access in the Adult Patient: 2019 Update. Peritoneal Dialysis International, 39(5), 414–436. https://doi.org/10.3747/pdi.2018.00232
- Litherland J, Lupton EW, Ackrill PA, et al. Computerized tomographic peritoneography: CT manifestations of in the investigation of leaks and abnormal collections in patients on CAPD. Nephrol Dial Transplant 1994;9:1449–52.

- Juergensen PH, Rizvi H, Claride VJ, et al. Value of scintigraphy in chronic peritoneal dialysis patients. Kidney Int 1999;55:1111–19.
- Prischl FC, Muhr T, Seiringer EM, et al. Magnetic resonance imaging of the peritoneal cavity among peritoneal dialysis patients, using the dialysate as "contrast medium". J Am Soc Nephrol 2002;13:197–203.

Baxter

